


Capsule Endoscopy Systems



Capsule Endoscopy Systems

Intended use

Diagnosis and screening of small bowel disease such as obscure gastrointestinal bleeding, Crohn's disease, celiac disease, Familial polyposis syndrome etc.

Features

- No waiting for image data importing.
- Real-time display monitor on receiving device.
- Less equipment, more efficient use.



Specifications

CapsuBot	CapsuRec	CapsuStation
 <p>Model: SX-CS Physical weight: ≤ 4.4g Material: polycarbonate Field of view: 160° Image resolution: 256x256 Frame rate: 2fps Battery duration: ≥ 13 hrs (Li-MnO₂) Voltage: 3 v Sealing: IPX8</p>	 <p>Model: SX-CSX-1 (data recorder) Model: SX-CSA-1 (antenna vest) Size: 135x71x23 mm Lithium-Ion Battery Voltage: 3.7v Duration: ≥ 14hrs Storage: SD Card (8 GB) Battery capacity: 5000mAh Weight: 310g \pm 10g Vest antennas: 8 pcs</p>	 <p>Compatible with: Windows XP, Windows 7, Windows 8, Windows 10; Easy installation Very efficient experience</p>